

CALCEIN STAINING

Protocol used for calcein bulk staining of microcracks

Created by Silvia Owusu (✉ silvia.owusu@zmk.unibe.ch) on September 29 2021
Created by Mathieu Simon (✉ mathieu.simon@artorg.unibe.ch) on September 29 2021

Last edit: September 23, 2022

Protocol purpose:

This protocol is aimed to stain in vivo microcracks by performing bulk staining of 3 mm thick bone samples using calcein. Calcein binds to free calcium present at the surface of microcracks and bulk staining allows to stain microcracks present in vivo only, not the ones generated by sample preparation.

PROTOCOL

Step 1 [~30 + 5 minutes minutes]

Cutting

- Cut bone sample using a diamond band saw into slices of 3 mm thickness.

Step 2 [~30 + 3 days minutes]

Fixation

- Put the samples in 4% formaldehyde for two days and rinsed for 8-12 hours in tap water. Exchanged water four times in the rinsing step.

Step 3 [~? + ? minutes]

Bulk staining

- Dissolve calcein at a concentration of 0.5 mM. Dissolution is performed either in distilled water or in 100% ethanol.
- Perform bulk staining under low vacuum (~300 mbar) at room temperature. 2 samples are stained for 7 hours and 2 are for 24 hours.

Step 4 [~? + ? minutes]

Embedding

- Do something

Step 5 [~? + ? minutes]

Sample preparation

- Do something

Step 6 [~? + ? minutes]

Histology analysis

- Do something

WARNINGS



Calcein

The $C_{30}H_{26}N_2O_{13}$ is bad for you. In case of contact, wash directly and inform the safety responsible.



Human tissue

Human bone can transmit diseases. Always use gloves when manipulating.

NOTES

REFERENCES